

LIFEPAK® 15 Monitor/Defibrillator: Pulse Oximetry and Pulse CO-Oximetry

SpO₂ / SpCO / SpMet

Connects from 15 directly to patient using reusable sensor



MASIMO® SET RAINBOW® DIRECT CONNECT SENSOR

3 ft Adult : 11996-000335
3 ft Pediatric: 11996-000337
12 ft Adult: 11996-000336
12 ft Pediatric: 11996-000338

*Adult for patients >30 kg
Pediatric for patients 10-50 kg*

SpO₂ Only

OPTION A
Connects from 15 directly to patient using reusable sensor



MASIMO SET RED™ LNCS DIRECT CONNECT SENSOR

3 ft Adult : 11996-000331
3 ft Pediatric: 11996-000333
12 ft Adult: 11996-000332
12 ft Pediatric: 11996-000334

*Adult for patients >30 kg
Pediatric for patients 10-50 kg*

OPTION B
Connects from 15 to patient using patient cable plus a reusable or disposable sensor



MASIMO SET LNCS REUSABLE SENSORS

LNCS DCIP Pediatric Digit Sensor
Reusable Sensor for Patients 10-50 kg
11171-000018

LNCS DCI Adult Digit Sensor
Reusable Sensor for Patients >30kg
11171-000017



MASIMO SET LNCS ADHESIVE DISPOSABLE SENSORS

LNCS Inf-L Infant Sensor
Single Patient Use for Patients
3-20 kg (20/box)
11171-000026

LNCS Neo-L Neonatal Sensor
Single Patient Use for Patients
<3 kg or >40 kg (20/box)
11171-000021

LNCS NeoPt-L Sensor
Single Patient Use for Patients
<1 kg (20/box)
11171-000022

LNCS Amtx Adult Sensor
Single Patient Use for Patients
> 30 kg (20/box)
11171-000019

LNCS Pmtx Pediatric Sensor
Single Patient Use for Patients
10-50 kg (20/box)
11171-000020



LIFEPAK® 12 Defibrillator/Monitor: Pulse Oximetry - SpO₂ Monitoring Only

MASIMO SET LNCS PATIENT CABLES

Connects from 12 to patient via patient cable plus a reusable or disposable sensor

MASIMO SET LNCS REUSABLE SENSORS

LNCS DCIP Pediatric Digit Sensor
Reusable Sensor for Patients 10-50 kg
11171-000018



LNCS DCI Adult Digit Sensor
Reusable Sensor for Patients >30 kg
11171-000017



MASIMO SET LNCS ADHESIVE DISPOSABLE SENSORS

LNCS Inf-L Infant Sensor
Single Patient Use for Patients
3-20 kg (20/box)
11171-000026



LNCS Neo-L Neonatal Sensor
Single Patient Use for Patients
<3 kg or >40 kg (20/box)
11171-000021



LNCS NeoPt-L Sensor
Single Patient Use for Patients
<1 kg (20/box)
11171-000022



LNCS Amtx Adult Sensor
Single Patient Use for Patients
> 30 kg (20/box)
11171-000019



LNCS Pmtx Pediatric Sensor
Single Patient Use for Patients
10-50 kg (20/box)
11171-000020



LNCS CABLES

4 ft 11171-000024
10 ft 11171-000016
14 ft 11171-000025



LNCS EXTENSION CABLE

4-foot extension patient cable (May be used with devices equipped with Nellcor conventional SpO₂)
11171-000027



For a complete listing of Nellcor® Sensors, LNOP® cables and MASIMO or Nellcor adapters, or for more information on Physio-Control products, visit store.physio-control.com or contact your sales representative.